

Tablet Container



Tablet container

for the storage and transport of tablets,
dragees and capsules.

GMP-compliant execution

Tablet containers are used for the storage and transport of pharmaceutical products. In line with the different customer-specific handling requirements to be met in regard to storing, the equipment and design of the containers is determined individually. A large range of accessories enables to meet with all demands. Different container uptakes, lifting lugs, transport rollers, paper bails and special stacking corners enable the simple and flexible handling of the containers at the place of operation. Covers and discharge valves can be determined individually.



MASCHINPEX
Maschinenbau GmbH

A Company in the
DRIAM-Group

Claude-Dornier-Straße 11
D-88239 Wangen im Allgäu
(Location: Industrial Park
Geiselharz-Schauwies)
Phone: (+49) 7520-9669-10
Fax: (+49) 7520-9669-11
Internet: www.maschinpex.de
E-Mail: info@maschinpex.de

*1000 l tablet container,
equipped with special
stacking corners for
space-saving storage.*



*Screwed butterfly
valve (Ø 200 mm).*



*Butterfly valve with
clamping ring for easy
dismounting and
cleaning of the
product outlet.*



200 l tablet container, equipped with lifting lugs for easy transport.

Technical Data Tablet-Container

Contents

Standard contents: 200, 250, 375, 500, 625 l,
or as desired

Materials

1.4404 or AISI 316L or
1.4301 or AISI 304

Surface

Ground, $R_a \leq 0,5 \mu\text{m}$,
or electropolished (on option)

Product outlet

Conical outlet with butterfly valve (Ø 150 mm),
or as desired

Product inlet

Manhole (Ø 400 mm) with clamping ring,
or as desired

Available on option

Electropolished surfaces
Guide rolls
Stacking corners
Container uptakes (lifting lugs, thorns, etc.)
Paper bails

The MASCHINPEX
product range includes:

- Inspection Machines
- Sorting Machines
- Outside Washing and
Decontamination Machines
- Capsule Polishers
- Feed Systems
- Containers